

Comparison of Fusion with VITRIOX GAS and VITRIOX ELECTRIC

▪ Introduction

Fusion based on gas and based on an electrical muffle oven differ considerably in the heating method. At the same time, the user expects comparable results.

▪ Technical Differences

In the case of the gas-based **VITRIOX GAS**, the crucible and the casting dishes are heated from below by a flame. The flame is formed from gas (propane, butane, natural gas) and oxygen (compressed air and oxygen). The relationship between gas/air/oxygen ultimately determines the temperature of the flame and therefore the temperature in the crucible. It is important that each burner is individually adjustable so that the temperature for each burner



Figure 1: VITRIOX GAS

can be set exactly. This is equally true for crucible burners as well as for die burners. There is no direct temperature measurement, but the gas flow is regulated so that even for a fluctuating inlet pressure, exactly the same amount of gas is always directed to the burner; leading to a higher temperature stability. The temperature is indirectly determined from the melting point of pure lithium tetraborate. The rule is that 9 g lithium tetraborate (FX-x100) melts in 4 minutes.

Comparison of Fusion with VITRIOX GAS and VITRIOX ELECTRIC



Figure 2: VITRIOX ELECTRIC

In the case of the **VITRIOX ELECTRIC** electrical fusion machine, the crucible and the casting dishes are evenly heated in a muffle oven. The muffle oven has a thermocouple that indicates the exact temperature. However, the temperature display must be factory calibrated, as the position of the thermocouple in relation to the position of the crucible is slightly different in each oven due to manufacturing tolerances. The goal is that every muffle oven displays the same temperature and also delivers the same analytical result.

▪ Illustration



Figure 4: VITRIOX GAS burners




Figure 3: VITRIOX ELECTRIC oven chamber

Comparison of Fusion with VITRIOX GAS and VITRIOX ELECTRIC

▪ Analytical Results

The following tables display the results for real customer samples that were fused with both the VITRIOX GAS gas fusion machine and the VITRIOX ELECTRIC electrical fusion machine. A separate calibration was used for each fusion system, which leads to slightly higher or lower means for one or another result.

**Table 1: Fusion of ferrochrome, analyzed with the
FLUXANA FERRO application package.**

Mass%	GAS	GAS	ELECTRIC	ELECTRIC
Preparation	No 1	No 2	No 1	No 2
Cr	65.75	65.55	66.16	65.63
Fe	23.95	23.33	23.18	23.98
Mn	0.17	0.16	0.11	0.14
Si	1.35	1.28	1.34	1.41

**Table 2: Fusion of slag, analyzed with the
FLUXANA SLAG application package.**

Mass%	GAS	GAS	ELECTRIC	ELECTRIC
Preparation	No 1	No 2	No 1	No 2
Al ₂ O ₃	10.97	10.99	10.79	10.79
CaO	0.64	0.65	0.63	0.63
Cr ₂ O ₃	4.44	4.45	4.60	4.63
Fe ₂ O ₃	0.99	0.99	0.99	0.99
K ₂ O	0.27	0.27	0.29	0.29
MgO	35.85	35.66	35.62	35.68
Mn ₂ O ₃	0.19	0.19	0.18	0.18
SiO ₂	46.23	46.12	46.04	46.04
SO ₃	0.17	0.17	0.19	0.18
SrO	0.14	0.14	0.12	0.12
TiO ₂	0.20	0.20	0.21	0.21

Comparison of Fusion with VITRIOX GAS and VITRIOX ELECTRIC

**Table 3: Fusion of quartzite, analyzed with the
FLUXANA RAW-PROF application package.**

Mass%	GAS	GAS	ELECTRIC	ELECTRIC
Preparation	No 1	No 2	No 1	No 2
Al ₂ O ₃	3.25	3.28	3.23	3.29
BaO	0.03	0.03	0.03	0.02
CaO	0.99	0.98	0.96	0.96
Cr ₂ O ₃	0.35	0.35	0.33	0.33
CuO	0.02	0.00	0.01	0.01
Fe ₂ O ₃	1.29	1.31	1.34	1.30
K ₂ O	0.98	0.98	1.00	1.02
MgO	0.26	0.25	0.28	0.30
Mn ₂ O ₃	0.04	0.04	0.03	0.03
Na ₂ O	0.55	0.52	0.42	0.46
P ₂ O ₅	0.07	0.06	0.07	0.07
SiO ₂	92.08	91.77	92.03	92.03
TiO ₂	0.22	0.23	0.23	0.23

**Table 4: Fusion of chrome ore, analyzed with the
FLUXANA RAW-PROF application package.**

Mass%	GAS	GAS	ELECTRIC	ELECTRIC
Preparation	No 1	No 2	No 1	No 2
Al ₂ O ₃	7.66	7.70	7.73	7.79
Cr ₂ O ₃	48.28	48.18	48.47	48.48
Fe ₂ O ₃	14.80	14.79	14.54	14.53
MgO	20.03	20.16	20.26	20.26
NiO	0.23	0.23	0.21	0.22
SiO ₂	7.81	7.81	7.70	7.72

Comparison of Fusion with VITRIOX GAS and VITRIOX ELECTRIC

Table 5: Fusion of cement, analyzed with the
FLUXANA CEM application package.

Mass%	GAS	GAS	ELECTRIC	ELECTRIC
Preparation	No 1	No 2	No 1	No 2
Al ₂ O ₃	5.58	5.58	5.54	5.55
CaO	64.63	64.68	64.77	64.77
Cr ₂ O ₃	0.01	0.01	0.01	0.01
Fe ₂ O ₃	2.61	2.61	2.66	2.65
K ₂ O	0.87	0.87	0.87	0.86
MgO	1.44	1.45	1.40	1.43
Mn ₂ O ₃	0.06	0.07	0.06	0.06
Na ₂ O	0.24	0.23	0.21	0.23
P ₂ O ₅	0.16	0.16	0.16	0.16
SiO ₂	20.92	20.96	20.83	20.86
SO ₃	2.95	2.91	3.05	3.07
SrO	0.28	0.28	0.29	0.29
TiO ₂	0.29	0.28	0.29	0.28
ZnO	0.05	0.05	0.05	0.05

▪ Summary

The analytical data shows that with optimized settings both the VITRIOX GAS gas fusion machine and the VITRIOX ELECTRIC electrical fusion machine provide excellent consistent results.

From experience, it can be said that the fusion time of both systems is also approximately the same.

Volatile elements such as the halogens F and Cl escape the fusion in the open system of the gas driven machine but can be determined with sufficient precision in the electrical fusion machine, due to its closed construction.

Another characteristic of the electrical system is the use of quartz crucibles instead of platinum crucibles. In this way, fusions for applications that strongly damage platinum crucibles (e.g., metals or sulfides) can also be performed.

Comparison of Fusion with VITRIOX GAS and VITRIOX ELECTRIC

Literature

- [1] Rainer Schramm, Röntgenfluoreszenzanalyse in der Praxis, korrigierte Auflage II, FLUXANA (2017).
- [2] FXBA_1014_01 White Paper VITRIOX ELECTRIC Chlorine in Cement
- [3] FXBA_1015_01_White_paper_Einsatz_von Quarztiegeln_als_Platinersatz
- [4] FXBA_1007_02_Applikationsbericht_Kalibrierung_von_Fluor_in_Zement_mittels_Boratschmelzen
- [2] www.fluxana.com